
53 64

pH

Key points. Hauptteile.
Puntos clave. Points clé.
Punti chiave. Pontos chaves.

- ① Protective cap.
Schutzkappe.
Tapón protector.
Bouchon protecteur.
Tappo protettore.
Tampa protetora.
- ② S8 connector.
S8-Stecker.
Conector S8.
Connecteur S8.
Connettore S8.
Conector S8.
- ③ PG 13.5 thread.
Montagegewinde PG 13,5.
Rosca de sujeción PG 13,5.
Filetage PG 13,5.
Filettatura PG 13,5.
Rosca de fixação PG 13,5.
- ④ Pressure compensator.
Druckkompensator.
Compensador de presión.
Compensateur de pression.
Compesatore di pressione.
Compensador de pressão.
- ⑤ Reference element.
Referenzelement.
Elemento de referencia.
Élément de référence.
Elemento di riferimento.
Elemento de referência.
- ⑥ Solid electrolyte.
Festelektrolyt.
Electrolito sólido.
Électrolyte solide.
Elettrolita solido.
Eletrólito sólido.
- ⑦ Diaphragm: open, 2 lateral.
Diaphragma: 2 seitliche Öffnungen.
Diafragma: abierto, 2 laterales.
Diaphragme: ouvert, 2 latérales.
Diaframma: aperto, 2 laterali.
Diafragma: aberto, 2 laterais.
- ⑧ Sensitive glass membrane.
Sensitive Glasmembran.
Membrana sensible.
Membrane sensible.
Membrana sensibile.
Membrana sensível.
- ⑨ Storage protector.
Lagerungsschutz.
Protector de almacenamiento.
Tube protecteur de stockage.
Protettore per la conservazione.
Protetor para armazenamento.

(EU Directive 2002/96/EC)

pH electrode 53 64

Introduction

The **53 64** is a low maintenance pH electrode for extremely difficult media.

SAFETY NOTE

Before use with hazardous substances, find out about the relevant hazards, protective measures, and what to do in an emergency.

Electrode preparation

- Remove the storage protector ⑨.
- Verify that there are not air bubbles inside the membrane ⑧. They can be removed by shaking the electrode as a clinical thermometer.

Installation

In immersion, insertion, flotation housings, etc. with PG 13,5 thread.

 See *Housing instruction manual*.

Cable. They are connected to coaxial cables with an AS9 connector.

Calibration

 See *pH-meter instruction manual*.

Measurements. Recommendations

See pH-meter instruction manual.

Ensure that the electrode's specifications are adapted to the process.

Sample must cover the diaphragm ⑦.

Maintenance

- Deposits forming on the electrodes, caused by continual exposure of the sensors to a dirty or contaminated medium, is one of the most frequent causes of problems when measuring pH. When possible it's recommended to use a housing with automatic sensor cleaning system.
- When an electrode gets contaminated, the simple rinsing with distilled water is not sufficient. The most suitable cleaning solution is one which acts in a more selective way over the contamination (look for regeneration solutions).
- **Re-hydration.** A "dry" electrode recovers by immersing it in diluted HCl solution for several hours. It is recommended to read the information supplied with the electrode's regenerating solutions.

Troubleshooting

... the meter always shows the same pH value.

Short circuit in the electrode or cable.

Substitute the cable or the electrode.

... shifting of pH values.

Reference system contamination.

Replace the electrode

... unstable reading.

Diaphragm blockage.

Clean, regenerate, replace the electrode.

... impossible to perform calibration.

... very slow response.

Dirt on the diaphragm or the membrane.

De-hydrated membrane.

Electrode ageing or wear.

Clean, regenerate, replace the electrode.

Storage

The connector should be protected from humidity by the protective cap ①.

Place the electrode in the storage protector ⑨ with KCl 3M.

Life of an electrode

The average life expectancy of a pH process electrode oscillates between 6 or 12 months, depending on the working conditions (in standard applications).

Electrodes exposed to high temperatures do not last as long. Very alkaline mediums also reduces the life span of electrodes.

 *If an electrode has too short life span, less than three months, ask us about it. Probably there is a more suitable electrode for the application.*

Specifications

Measuring range	0 to 14 pH
Operating temp.	0 to 130 °C
Max. pressure	6 bar
Reference element	cartridge
Diaphragm	open, 2 lateral
Electrolyte	solid
Body material	glass
Thread	PG 13,5
Connector	S8
Min. immersion depth	15 mm 

Applications

Samples with dyes, pigments and colourings. Fermentation. Sugar industry. Samples with oil. Organic solvents.

Resistant to sterilisation with steam and in autoclave.

Resistant to CIP processes

CERTIFICATE OF QUALITY

The attached electrode 53 64 has passed the quality test and follows the next specifications:

- Asymmetry potential < ± 15 mV.
- Sensitivity, pH 4 to 7 (at 25 °C) > 98%.
- Response time, pH 4 to 7 < 20 s.

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