



# Biotrate Digital Burette

Convenient Titration

Simplifying Progress

**SARTORIUS**



# Biotrate Premium Digital Burette

## Convenient Titration

Sartorius Biotrate is a premium digital burette equipped with sophisticated functions for easy titration. Biotrate offers high chemical resistance and it is a smart solution for titrating various liquids accurately and safely. Biotrate is easy and convenient to use due to its large and clear electronic display, smooth operating wheels, and 360° rotating dispensing head.

### High Chemical Resistance

- Highly resistant parts in the liquid pathway guarantee excellent chemical resistance
- The softly-moving piston ensures reliable results

### Premium Quality

- High quality parts and materials guarantee highly accurate and precise results as well as a long lifetime

### Convenient to Use

- Biotrate is comfortable to use due to soft touch operating wheels and 360° rotating dispensing head
- Extremely smooth and light to operate

## Biotrate Features

### Chemical Resistance

Biotrate offers excellent chemical resistance due to its high-quality materials, allowing the use of a wide variety of liquids.

### Versatile Accessories

Biotrate can be equipped with a wide range of accessories, e.g., various bottle adapters, dispensing heads, and aspiration tubes.

### Light to Use

The operating wheels are light and comfortable to use due to the soft touch surface and smoothly moving piston. These guarantee good ergonomics and excellent results in everyday use.

### Long Battery Life

Biotrate is powered with a long-lasting battery so there is no need to plug in or recharge the unit. The battery is also easy to replace by the user.

### Large and Clear Display

The display clearly indicates the measured volume, and shows whether you are measuring an aspirated or titrated volume.

### Easy to Clean

Biotrate is easy to disassemble for cleaning and autoclaving. The liquid pathway is fully autoclavable.

### User Adjustment

Biotrate can be adjusted for various liquids and it is easy to revert to factory settings. The display also indicates when a custom adjustment is in use.

### Find out more

For more information, please visit  
[www.sartorius.com/en/products/pipetting/bottle-top-dispensers](http://www.sartorius.com/en/products/pipetting/bottle-top-dispensers)

# Convenient to Use



## Smart 360° Rotating Dispensing Head

The dispensing head can be turned to any position without having to turn the bottle or the entire burette. The display, operating wheels, and the bottle label remain in optimal position ensuring comfortable use.



## Media Recirculation System

Biotrate is equipped with a media recirculation system that allows unused media to be recovered from the device and returned to the bottle. This ensures easy priming and minimizes the media loss.



## Connector for Drying Tube and Filters

Drying tubes and filters can be connected to the back of the Biotrate. It is recommended that drying tubes be used when titrating humidity-sensitive liquids.



## UV-Protection Window

Biotrate is supplied with an amber coloured replacement window for protecting light-sensitive media.



## Easy to Attach Aspiration Tube

The aspiration tube is easy to attach with a secure threaded connection. This prevents air from entering the system and ensures reliable results.



## Easily Replaceable Valves

Valves can be easily replaced by the user so there is no need to send the unit to be serviced.

# Sartorius Minisart® Syringe Filters

Sartorius Minisart® filters fit Biotrate digital burettes perfectly. Minisart® filters offer superior flow rate, high total throughput volumes, and exceptional quality.



Biotrate can be equipped with a Minisart® filter for media filtration. Dispensing heads with luer lock connection are available separately.



For filtering replacement air, Biotrate can be equipped with a Minisart® filter.

# Sartorius Pipettes, Pipette Tips, and Services

Our leading range of superior electronic and mechanical pipettes are the perfect solution for your liquid handling needs. Sartorius pipettes and tips fit seamlessly together, maximizing pipetting performance and ensuring accuracy and precision. For more information, visit [www.sartorius.com](http://www.sartorius.com).

Sartorius also provides pipette maintenance, repair, and calibration services globally. Contact your nearest Sartorius service center or office for more information.



# Ordering Information

## Biotrate

### Prospenser

Order Code	Item	Nominal Volume
LH-723080	Biotrate	10 mL
LH-723081	Biotrate	20 mL
LH-723082	Biotrate	50 mL

Note: All models have A45 thread as a default and they are supplied with A32, A38, and S40 bottle adapters. Bottle is not included.

### Accessories

Order Code	Item	Material
LH-721650	Dispensing head, standard	FEP
LH-721651	Dispensing head, fine tip	FEP
LH-721648	Dispensing head, 1.5 mm spiral	FEP
LH-721652	Dispensing head, 0.4 m adjustable	FEP
LH-721653	Dispensing head, luer lock	FEP/PP
LH-721654	Dispensing head, luer lock	FEP/PFA
LH-721678	Aspiration tube, 310 mm	FEP
LH-721679	Aspiration tube, 400 mm	FEP
LH-721682	Recirculation tube, 70 mm	FEP
LH-721689	Bottle adapter, A28	ETFE
LH-721684	Bottle adapter, A28	PP
LH-721683	Bottle adapter, A32	PP
LH-721688	Bottle adapter, A32	ETFE
LH-721686	Bottle adapter, A38	PP
LH-721733	Bottle adapter, A38	ETFE
LH-721685	Bottle adapter, A38/430	PP
LH-721732	Bottle adapter, A38/430	ETFE
LH-721687	Bottle adapter, S40	PP
LH-721734	Bottle adapter, S40	ETFE
LH-721735	Bottle adapter, NS29/32	Silicone
LH-721743	Bottle support, for bottles 75–120 mm	Silicone



ISO9001  
13485  
14001  
Sartorius  
Certified






**Germany**

Sartorius Lab Instruments GmbH & Co. KG  
Otto-Brenner-Strasse 20  
37079 Goettingen  
Phone +49 551 308 0

**USA**

Sartorius Corporation  
565 Johnson Avenue  
Bohemia, NY 11716  
Phone +1 631 254 4249  
Toll-free +1 800 635 2906

 For additional information,  
visit [www.sartorius.com](http://www.sartorius.com)