

Corning Incorporated  
Life Sciences

Registered  
ISO 9001:2008

## Product Description

**Catalog Number:** 1300P-1L

**Product Description:** Corning® Bottle, Washing, Reusable Plastic, LDPE, Not Sterile, 1000 mL  
GL-32 PP Cap and PP Spray Tube

**Component Materials:**

Bottle - Low Density Polyethylene (LDPE)  
Stopper - Polypropylene (PP)  
Spray Tube - Polypropylene (PP)

**Attributes:**

Capacity - 1000 mL  
Product Dimensions - 92 mm Body OD x 221 mm High without Cap

**Autoclaving:**

This product IS NOT suitable for Autoclaving.

**Heating / Cooling:**

The upper temperature limit for this plastic material is 80°C.

This product IS NOT suitable for Dry Heat Sterilization as it will melt. This product MUST NOT be heated by any direct, surface contacting heat source.

The lower limit usage temperature for this product is -50°C.

**Vacuum / Pressure:**

This product is not designed for and is not rated for vacuum or pressure applications.

**Sterilization:**

Shipped Not Sterile

**Lot Number Designation:**

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

**Revision Date:**

01 OCT 2012

Rev No: 1